



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Brand UV cuvette micro, center height 8.5 mm,  
filling volume 70 µl to 850 µl**

**€132.75** As low  
as  
**€119.75**  
plus VAT &  
Shipping

## Product Images



## Description

---

Specially developed for the photometric determination of proteins, ssDNA, dsDNA, RNA and oligonucleotides in the UV range. Applicable for wavelengths from 230 nm, ideally suited for measurements at 260 nm, 280 nm and in the visible wavelength range. Standard layer thickness 10 mm. Even small sample quantities from 70 µl are sufficient for measurement. Individually packaged UV cuvettes micro are DNase, DNA and RNase free

## Additional Information

---

No.	BN-2101
Manufacturer (Brand)	Brand
VGKL number	294959210
EAN	4058072213282
Gross weight	1.1951
HSNumber	90279000
Transport temperature	Room temperature
Storage temperature	Room temperature
Net weight	0.75350 kg
Country of origin	Germany
Packing width	0.22700 m
Packing height	0.20500 m
Packing depth	0.36000 m
Packaging volume	16752.60000000 ccm
Customs tariff number	90279000

