



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Brand UV cuvettes semi-micro, 1.5-3.0 ml, 100 pcs./pack

€38.90
plus VAT &
Shipping

Product Images



Description

made of plastic, high chemical resistance, nest number sorted.

Wavelength range: 220–900 nm

Standard deviation in absorbance units:

240 nm $\leq \pm$ 0.007

300 nm $\leq \pm$ 0.005

Layer thickness: 10 mm

UV cuvette semi-micro

Filling volume: 1.5 to 3.0 ml

Window dimensions: 4.5 x 23 mm (min.)

Packaging unit: 100 pcs./pack

Additional Information

No.	BN-2098
Manufacturer (Brand)	Brand
old neoLab article no.	1-6250
EAN	4058072213251
Transport temperature	Room temperature
Material	Polystyrene (PS)
Volume MAX	3 ml
Volume MIN	1.5 ml
Cuvette properties	with special chemical resistance for 1xuse suitable for UV range
Layer thickness	10.0 mm
Standard deviation	at 240 nm $< \pm$ 0.007 at 300 nm $< \pm$ 0.005
Type cuvette	Semi-micro
Wavelength	220 to 900 nm

