



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Brand UV cuvettes micro 8.5 mm, 70-850 µl, 100 pcs/pack

€29.15 As low as
€26.55
**plus VAT &
Shipping**

Product Images



Description

made of plastic, high chemical resistance, nest number sorted.

Wavelength range: 220–900 nm

Standard deviation in absorbance units:

240 nm $\leq \pm$ 0.007

300 nm $\leq \pm$ 0.005

Layer thickness: 10 mm

UV cuvette micro: from 70 μ l

Round lid for secure closure and storage down to -20°C .

Individually packed UV cuvettes micro are DNase-, DNA- and RNase-free!

Filling volume: 70 to 850 μ l

Dimensions window: 2 x 3.5 mm (min.)

Packaging unit: 100 pcs./pack

Additional Information

No.	BN-2100
Manufacturer (Brand)	Brand
old neoLab article no.	1-6221
EAN	4058072213275
Transport temperature	Room temperature
Material	Polystyrene (PS)
Volume MAX	0.85 ml
Volume MIN	0.07 ml
TBST MIN	-20°C
Cuvette properties	for 1x use suitable for UV range cap/lid available as accessory
Layer thickness	10.0 mm
Wavelength	220 to 900 nm
Standard deviation	at 240 nm $< \pm$ 0.007 at 300 nm $< \pm$ 0.005
Type cuvette	Micro
Center height	8.5 mm

