



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Brand UV cuvettes micro 15 mm, 70-550 µl,  
individually packed, 100 pcs./pack**

**€60.55  
plus VAT &  
Shipping**

## Product Images



## Description

---

made of plastic, high chemical resistance, nest number sorted.

Wavelength range: 220–900 nm

Standard deviation in absorbance units:

240 nm  $\leq \pm$  0.007

300 nm  $\leq \pm$  0.005

Layer thickness: 10 mm

UV cuvette micro: from 70  $\mu$ l

Round lid for secure closure and storage down to  $-20^{\circ}\text{C}$ .

Individually packed UV cuvettes micro are DNase-, DNA- and RNase-free!

Filling volume: 70 to 550  $\mu$ l

Dimensions window: 2 x 3.5 mm (min.)

Packaging unit: 100 pcs./pack, individually packed

## Additional Information

---

No.	BN-2105
Manufacturer (Brand)	Brand
old neoLab article no.	1-6219
EAN	4058072213329
Transport temperature	Room temperature
Material	Polystyrene (PS)
Volume MAX	0.5 ml
Volume MIN	0.07 ml
TBST MIN	$-20^{\circ}\text{C}$
Cuvette properties	for 1x use   suitable for UV range   individually packed   cap/lid available as accessory   RNase-free   DNase-free   DNA-free
Layer thickness	10.0 mm
Standard deviation	at 240 nm $< \pm$ 0.007   at 300 nm $< \pm$ 0.005
Type cuvette	Micro
Wavelength	220 to 900 nm

