



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Hellma Ultra Micro Cell 105.201-QS SD 10mm, VOL 100µl, Z.15mm

**€376.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Ultra Micro Cell Type 105.201-QS with stopper made of PE for absorbance measurements, window material High Performance Quartz Glass (200 nm - 2500 nm), volume 100 µl, path length 10 mm, center height 15 mm, outer dimensions 45 x 12.5 x 12.5 mm, aperture 5 x 2 mm

## Additional Information

---

No.	KI-0488
Reg. trademark	SUPRASIL®
Manufacturer (Brand)	Hellma
Dual use indicator	false
VGKL number	105512011
EAN	4058072785444
Gross weight	0.03
Transport temperature	Room temperature
Net weight	0.01150 kg
Material	Quartz glass
Volume MAX	0.1 ml
Country of origin	Germany
Region of origin	Rhineland-Palatinate
Packing width	0.04000 m
Packing height	0.02500 m
Packing depth	0.06400 m
Customs tariff number	90279050
Cuvette properties	for spectroscopy
Layer thickness	10.0 mm
Wavelength	200 to 2500 nm
Center height	15 mm

