



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Hellma Ultra Micro Cell 105.250-QS, SD 10x2mm, VOL 100µl, Z.15mm

**€501.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Ultra-micro cuvette type 105.250-QS with stopper made of PE for fluorescence measurements, window material High Performance Quartz glass (200 nm - 2500 nm), number of windows 3, volume 100 µl, path length 10 x 2 mm, center height 15 mm, outer dimensions 45 x 12.5 x 12.5 mm, aperture 5 x 2 mm

## Additional Information

---

No.	KI-0495
Reg. trademark	SUPRASIL®
Manufacturer (Brand)	Hellma
VGKL number	105512501
Dual use indicator	false
EAN	4058072785505
Gross weight	0.03
Transport temperature	Room temperature
Net weight	0.01120 kg
Material	Quartz glass
Volume MAX	0.1 ml
Country of origin	Germany
Region of origin	Rhineland-Palatinate
Packing width	0.04000 m
Packing height	0.02500 m
Packing depth	0.06400 m
Customs tariff number	90279050
Wavelength	200 to 2500 nm
Center height	15 mm

