



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Hellma Cylinder Macro Cell 120-QS SD 1mm, VOL 280µl

€261.00
plus VAT &
Shipping

Product Images



Description

Cylindrical cell type 120-QS with stopper made of PTFE for absorbance measurements, window material High Performance Quartz glass (200 nm - 2500 nm), volume 280 µl, path length 1 mm, outer diameter 22 mm, inner diameter 19 mm, outer depth 3.5 mm

Additional Information

No.	KI-0571
Reg. trademark	SUPRASIL®
Dual use indicator	false
Manufacturer (Brand)	Hellma
VGKL number	195516002
EAN	4058072786090
Gross weight	0.02
Transport temperature	Room temperature
Net weight	0.00350 kg
Material	Quartz glass
Volume MAX	0.28 ml
Country of origin	Germany
Region of origin	Rhineland-Palatinate
Packing width	0.04000 m
Packing height	0.02500 m
Packing depth	0.06400 m
Customs tariff number	90279050
Cuvette properties	for spectroscopy
Layer thickness	1.0 mm
Wavelength	200 to 2500 nm

