



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Hellma Cylinder Macro Cell 120-QS SD 2mm, VOL 560µl

€259.00
plus VAT &
Shipping

Product Images



Description

Cylindrical cell type 120-QS with stopper made of PTFE for absorbance measurements, window material High Performance Quartz glass (200 nm - 2500 nm), volume 560 µl, path length 2 mm, outer diameter 22 mm, inner diameter 19 mm, outer depth 4.5 mm

Additional Information

No.	KI-0572
Reg. trademark	SUPRASIL®
Manufacturer (Brand)	Hellma
VGKL number	195516012
EAN	4058072786106
Gross weight	0.02
Transport temperature	Room temperature
Material	Quartz glass
Volume MAX	0.56 ml
Packing width	0.04000 m
Packing height	0.02500 m
Packing depth	0.06400 m
Customs tariff number	90279050
Cuvette properties	for spectroscopy
Layer thickness	2.0 mm
Wavelength	200 to 2500 nm

