



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Hellma Cylindrical Cell 540.135-QS, VOL 14000 μ l

€205.00
plus VAT &
Shipping

Product Images



Description

Cylindrical cell type 540.135-QS with PTFE stopper for scattered light measurements, window material High Performance Quartz glass (200 nm - 2500 nm), volume 14000 µl, outer dimensions 75 x 20 mm, inner dimensions 74 x 18 mm

Additional Information

No.	KI-0658
Reg. trademark	SUPRASIL®
Manufacturer (Brand)	Hellma
VGKL number	105515405
Dual use indicator	false
EAN	4058072786892
Gross weight	0.02
Transport temperature	Room temperature
Net weight	0.01380 kg
Material	Quartz glass
Volume MAX	14 ml
Country of origin	Germany
Region of origin	Rhineland-Palatinate
Packing width	0.02000 m
Packing height	0.07500 m
Packing depth	0.02000 m
Customs tariff number	90279050
Cuvette properties	for spectroscopy
Layer thickness	20.0 mm
Wavelength	320 to 2500 nm

