



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Hellma glass filter set 666-S200

€1,150.00
plus VAT &
Shipping

Product Images

Description

Glass filter set type 666-S200 for checking photometric accuracy and linearity, consisting of 3 neutral glass filters (F2, F4, F203) with absorbances: 0.25 / 1.0 / 2.0, certified at wavelengths: 440; 465; 546.1; 590; 635 nm, in high-quality storage box with DAkKS calibration certificate Attention: Article available until 30.06.2022. Further information will follow.

Additional Information

No.	KI-0762
Manufacturer (Brand)	Hellma
EAN	4058072787271
Transport temperature	Room temperature
Country of origin	Germany
Customs tariff number	90279050

