




neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® Sterile Disposable CuveTips™ 6 x racks per carton

**€758.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Disposable reflector tips, the so-called CuveTip®, are used for measurement with the Photopette®. In contrast to conventional photometer measurements, only a small volume of the analyte is required. The liquid to be measured flows into the CuveTips® by itself; active sample transfer into the cuvette is not necessary. A minimum volume of 40 µl is required to fill the cavity of the CuveTip®. This volume can be recovered after the measurement.

## Additional Information

---

No.	C-3208
Manufacturer (Brand)	neoLab
EAN	4058072160319
Transport temperature	Room temperature

