




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Sterile Disposable CuveTips™ 12 x racks per carton

€1,247.00
plus VAT &
Shipping

Product Images



Description

Disposable reflector tips, the so-called CuveTip®, are used for measurement with the Photopette®. In contrast to conventional photometer measurements, only a small volume of the analyte is required. The liquid to be measured flows into the CuveTips® by itself; active sample transfer into the cuvette is not necessary. A minimum volume of 40 µl is required to fill the cavity of the CuveTip®. This volume can be recovered after the measurement.

Additional Information

No.	C-3210
Manufacturer (Brand)	neoLab
EAN	4058072160333
Transport temperature	Room temperature

