



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Eppendorf Deepwell Plate 96/1000 $\mu$ L, wells clear, 1,000 $\mu$ L, PCR clean, white, 20 plates (5 bags $\times$ 4 plates)

Special Price

**€129.00** was

**€172.00**

**plus VAT &  
Shipping**

### Product Images



## Description

---

Eppendorf Deepwell Plates 96 and 384 are high-performance plates for all manual and automated applications - from sample storage at -86 °C to DNA denaturation at 100 °C. Made of pure polypropylene (PP), the plates offer innovative features: More reliability. More efficiency. More stability. Plates are available with a white, yellow, green or blue rim. High-contrast alphanumeric labeling is laser-engraved into the rim for easy identification of wells. Find your sample up to 30% faster and reduce pipetting errors with this OptiTrack® matrix. With all their thoughtful features, these plates set new standards in the lab.

## Additional Information

No.	VB-0410
Manufacturer (Brand)	Eppendorf
VGKL number	106901209
EAN	4058072173005
Gross weight	2.3
HSNumber	39269097
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	2.07000 kg
Color	White
Material	Polypropylene (PP)
Volume	1000 µL
Country of origin	Germany
Packing width	0.64000 m
Packing height	0.26000 m
Wide	127.8 mm
Packing depth	0.46000 m
Packaging volume	77000.000000 ccm
Customs tariff number	39269097
Height	44.1 mm
Depth	85.5 mm
Number of cavities	96
Degree of purity	PCR clean
Shape floor	V-shape (conical)
centrifugable up to	6000 g

