



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Eppendorf Deepwell plate 96/500 $\mu$ l, Normal Pack, DNA LoBind, white, 40 plates (5 bags of 8)

Special Price  
**€312.00** was  
**€416.00**  
plus VAT &  
Shipping

### Product Images



## Description

---

DNA LoBind Plates improve the recovery of nucleic acids by reducing their binding to the vessel wall. A combination of special manufacturing technologies and selected polypropylene batches ensures nearly 100% recovery of DNA/RNA molecules - without surface coating to eliminate sample contamination hazards. DNA LoBind Plates are ideal for sample preparation or long-term storage of nucleic acids in forensic applications, microarrays and next-generation sequencing.

## Additional Information

No.	VB-0430
Manufacturer (Brand)	Eppendorf
VGKL number	106903104
EAN	4058072168445
Gross weight	2.75
HSNumber	39269097
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	2.47500 kg
Color	White
Material	Polypropylene (PP)
Volume	500 µL
Country of origin	Germany
Packing width	0.64000 m
Packing height	0.24000 m
Wide	127.8 mm
Packing depth	0.46000 m
Packaging volume	71000.000000 ccm
Customs tariff number	39269097
Height	27.1 mm
Depth	85.5 mm
Number of cavities	96
Degree of purity	PCR clean

