



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Brand strip pl.w.grid 96-well (12 str.xF8), PS,  
transp. 350 µl, immunoGrade HB; Packing unit:  
100 pieces (4 x 25 pcs.)**

**€655.60  
plus VAT &  
Shipping**

## Product Images



## Description

---

**BRANDplates® micro** plates have ideal optical properties and are manufactured according to ANSI/SLAS standards. They are also available with immunoGrade™ surface as 96-well **strip plates** - for high flexibility and optimal results.

- immunoGrade™ high binding capacity (highbinding), ideal for IgG and molecules with hydrophobic and hydrophilic domains

**They offer high flexibility. Only as many wells are used as needed.**

### Frame with grid:

12 strips of 8 wells each. Individual strips can be removed from the plate and reused. Strips can be separated into up to 8 individual wells and inserted individually into the frame plate.

## Additional Information

---

No.	BN-2386
Manufacturer (Brand)	Brand
old neoLab article no.	7-5392
EAN	4058072215934
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	white transparent (crystal clear)
Material	Polystyrene (PS)
Volume	350 µl
sterile	No
Multi-hard disk features	surface optimized alphanumeric coded
Number of cavities	96
Fluid behavior	amphiphil
suitable for cooling in	Refrigerator freezer/ freezer
Shape floor	F shape (flat)
Shape cavities	around
Type plate	Microtiter plate
Type frame MTP	full frame

