



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e


Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Brand microtiter plate, immunoGrade, 384-well,  
PS, standard, transp., F-B., 100 µl, VE=50 (10 x 5  
pcs.)**

**€273.75** As low  
as  
**€259.95**  
plus VAT &  
Shipping

## Product Images



## Description

---

For optimal binding of IgG.

- Optimized for binding IgG, highest binding capacity for molecules with hydrophilic and hydrophobic regions.
- The surface of choice for the majority of standard ELISA.
- Suitable for solid phase immunoassays.
- Analogous to high-binding products from other manufacturers.
- Free of endotoxins, DNase, DNA, RNase, non-cytotoxic.

## Additional Information

---

|                          |   |
|--------------------------|---|
| No.                      | BN-2335   |
| Manufacturer (Brand)     | Brand   |
| old neoLab article no.   | 7-5296  |
| EAN                      | 4058072215422                                     |
| Transport temperature    | Room temperature                                  |
| autoclavable at 121°C    | No  |
| Color                    | transparent (crystal clear)                       |
| Material                 | Polystyrene (PS)                                  |
| Volume                   | 1000 µL   |
| sterile                  | No  |
| Multi-hard disk features | alphanumeric coded RNase-free DNase-free DNA-free |
| Number of cavities       | 384   |
| suitable for cooling in  | Refrigerator freezer/ freezer                     |
| Shape floor              | F shape (flat)                                    |
| Shape cavities           | angular   |
| Type plate               | Microtiter plate                                  |
| Type frame MTP           | without frame                                     |

