



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Brand microtiter plate, immunoGrade, 384-well,
PS,standard, white, F-B., 100 µl, VE=50 (10 x 5
pcs.)**

€287.70 As low
as
€273.05
**plus VAT &
Shipping**

Product Images



Description

For optimal binding of IgG.

- Optimized for binding IgG, highest binding capacity for molecules with hydrophilic and hydrophobic regions.
- The surface of choice for the majority of standard ELISA.
- Suitable for solid phase immunoassays.
- Analogous to high-binding products from other manufacturers.
- Free of endotoxins, DNase, DNA, RNase, non-cytotoxic.

Additional Information

No.	BN-2336
Manufacturer (Brand)	Brand
old neoLab article no.	7-5297
EAN	4058072215439
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	White
Material	Polystyrene (PS)
Volume	1000 µL
sterile	No
Multi-hard disk features	alphanumeric coded RNase-free DNase-free DNA-free
Number of cavities	384
suitable for cooling in	Refrigerator freezer/ freezer
Shape floor	F shape (flat)
Shape cavities	angular
Type plate	Microtiter plate
Type frame MTP	without frame

