



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e


Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Brand microtiter plate, immunoGrade, 96-well,  
PS, standard, transp., F-B., 350 µl, PU=100 (20 x 5  
pcs.)**

**€166.95** As low  
as  
**€158.50**  
plus VAT &  
Shipping

## Product Images



## Description

---

For optimal binding of IgG.

- Optimized for binding IgG, highest binding capacity for molecules with hydrophilic and hydrophobic regions.
- The surface of choice for the majority of standard ELISA.
- Suitable for solid phase immunoassays.
- Analogous to high-binding products from other manufacturers.
- Clearly distinguishable by color code: blue, embossed alphanumeric coding on 96-well standard plates.
- Free of endotoxins, DNase, DNA, RNase, non-cytotoxic.

## Additional Information

---

No.	BN-2330
Manufacturer (Brand)	Brand
old neoLab article no.	7-5288
EAN	4058072215378
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	transparent (crystal clear)
Material	Polystyrene (PS)
Volume	350 µl
sterile	No
Multi-hard disk features	alphanumeric coded RNase-free DNase-free DNA-free
Number of cavities	96
suitable for cooling in	Refrigerator freezer/ freezer
Shape floor	F shape (flat)
Shape cavities	around
Type plate	Microtiter plate
Type frame MTP	without frame

