



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Brand Microtiterpl., cellGrade, 384-well, PS,
Standard, transp., F-B., 100 µl, VE=50
(individually packed, with lid)**

€411.10 As low
as
€390.50
plus VAT &
Shipping

Product Images



Description

For the cultivation of adherent cell cultures.

- Standard plate for the cultivation of adherent cell cultures.
- PS surface with various freely accessible chemical groups, such as carboxyl and hydroxyl groups.
- Surface is hydrophilic compared to untreated PS.
- Serum components are bound to the free chemical groups allowing indirect adhesion of cells.
- Free of endotoxins, DNase, DNA, RNase, non-cytotoxic.

Additional Information

No.	BN-2358
Manufacturer (Brand)	Brand
old neoLab article no.	7-5350
EAN	4058072215651
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	transparent (crystal clear)
Material	Polystyrene (PS)
Volume	100 µl
sterile	Yes
Multi-hard disk features	individually packed alphanumerically coded RNase-free DNase-free DNA-free
Number of cavities	384
Shape floor	F shape (flat)
Shape cavities	angular
Type plate	Microtiter plate
Type frame MTP	without frame
suitable for cooling in	Refrigerator freezer/ freezer

