



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e


Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Brand Microtiterpl., cellGrade, 384-well, PS,  
standard, white, F-B., 100 µl, VE=50 (individually  
packed, with lid)**

**€430.45** As low  
as  
**€408.95**  
plus VAT &  
Shipping

## Product Images



## Description

---

For the cultivation of adherent cell cultures.

- Standard plate for the cultivation of adherent cell cultures.
- PS surface with various freely accessible chemical groups, such as carboxyl and hydroxyl groups.
- Surface is hydrophilic compared to untreated PS.
- Serum components are bound to the free chemical groups allowing indirect adhesion of cells.
- Free of endotoxins, DNase, DNA, RNase, non-cytotoxic.

## Additional Information

---

No.	BN-2359
Manufacturer (Brand)	Brand
old neoLab article no.	7-5351
EAN	4058072215668
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	White
Material	Polystyrene (PS)
Volume	100 µl
sterile	Yes
Multi-hard disk features	individually packed alphanumerically coded RNase-free DNase-free DNA-free
Number of cavities	384
Shape floor	F shape (flat)
Shape cavities	angular
Type plate	Microtiter plate
Type frame MTP	without frame
suitable for cooling in	Refrigerator freezer/ freezer

