



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Brand microtiter plate, cellGrade, 96-well, PS,
transp.B., white, F-B., 330 µl, PU = 5 pcs**

**€66.75
plus VAT &
Shipping**

Product Images



Description

For the cultivation of adherent cell cultures.

- Standard plate for the cultivation of adherent cell cultures.
- PS surface with various freely accessible chemical groups, such as carboxyl and hydroxyl groups.
- Surface is hydrophilic compared to untreated PS.
- Serum components are bound to the free chemical groups allowing indirect adhesion of cells.
- Clearly distinguishable by color code: orange, embossed alphanumeric coding on 96-well standard plates.
- Free of endotoxins, DNase, DNA, RNase, non-cytotoxic.

Additional Information

No.	BN-2356
Manufacturer (Brand)	Brand
old neoLab article no.	7-5348
EAN	4058072215637
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	White
Material	Polystyrene (PS)
Volume	330 µl
sterile	Yes
Multi-hard disk features	individually packed alphanumerically coded RNase-free DNase-free DNA-free
Number of cavities	96
suitable for cooling in	Refrigerator freezer/ freezer
Shape floor	F shape (flat)
Shape cavities	around
Type plate	Microtiter plate
Type frame MTP	without frame

