



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Brand microtiter plate, cellGrade, 96-well, PS,  
transp.B., black, F-B., 330 µl, PU=5 pcs**

**€66.75  
plus VAT &  
Shipping**

## Product Images



## Description

---

For the cultivation of adherent cell cultures.

- Standard plate for the cultivation of adherent cell cultures.
- PS surface with various freely accessible chemical groups, such as carboxyl and hydroxyl groups.
- Surface is hydrophilic compared to untreated PS.
- Serum components are bound to the free chemical groups allowing indirect adhesion of cells.
- Clearly distinguishable by color code: orange, embossed alphanumeric coding on 96-well standard plates.
- Free of endotoxins, DNase, DNA, RNase, non-cytotoxic.

## Additional Information

---

No.	BN-2357
Manufacturer (Brand)	Brand
old neoLab article no.	7-5349
EAN	4058072215644
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	Black
Material	Polystyrene (PS)
Volume	330 µl
sterile	Yes
Multi-hard disk features	individually packed alphanumerically coded RNase-free DNase-free DNA-free
Number of cavities	96
Shape floor	F shape (flat)
Shape cavities	around
Type plate	Microtiter plate
Type frame MTP	without frame
suitable for cooling in	Refrigerator freezer/ freezer

