



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Brand microtiter plate, hydroGrade, 96-well, PS,
standard, transp., F-B., 350 µl, PU=100 (20 x 5
pcs.)**

€200.20 As low
as
€190.10
plus VAT &
Shipping

Product Images



Description

For binding hydrophilic molecules.

- Strongly hydrophilic, high affinity to hydrophilic molecules, such as glycoproteins and peptides, antibodies with predominantly hydrophilic regions and nucleic acids.
- Interesting alternative to immunoGrade™ surface for solid phase assays.
- Alternative for homogeneous assays with hydrophobic molecules, which then remain in solution.
- Clearly distinguishable by color code: blue, embossed alphanumeric coding on 96-well standard plates.
- Free of endotoxins, DNase, DNA, RNase, non-cytotoxic.

Additional Information

No.	BN-2339
Manufacturer (Brand)	Brand
old neoLab article no.	7-5300
EAN	4058072215460
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	transparent (crystal clear)
Material	Polystyrene (PS)
Volume	350 µl
sterile	No
Multi-hard disk features	alphanumeric coded RNase-free DNase-free DNA-free
Number of cavities	96
suitable for cooling in	Refrigerator freezer/ freezer
Shape floor	F shape (flat)
Shape cavities	around
Type plate	Microtiter plate
Type frame MTP	without frame

