



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Brand microtiter plate, inertGrade, 96-well, PS,  
transp.B., black, F-B., 330 µl, VE=40 (individually  
packed, with lid)**

**€1,015.30** As  
low as  
**€966.95**  
plus VAT &  
Shipping

## Product Images



## Description

---

For the cultivation of suspension cell cultures

- Particularly suitable for cell cultures when cell adhesion is not desired.
- Optimized surface properties reduce cell attachment, protein absorption, enzyme activation and cellular activation to a minimum.
- Stem cells can be prevented from premature differentiation, for example.
- Clearly distinguishable by color coding: orange, embossed alphanumeric coding on 96-well standard plates.
- Free of endotoxins, DNase, DNA, RNase, non-cytotoxic.

## Additional Information

---

No.	BN-2345
Manufacturer (Brand)	Brand
old neoLab article no.	7-5337
EAN	4058072215521
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	Black
Material	Polystyrene (PS)
Volume	330 µl
sterile	Yes
Multi-hard disk features	individually packed alphanumerically coded RNase-free DNase-free DNA-free
Number of cavities	96
suitable for cooling in	Refrigerator freezer/ freezer
Shape floor	F shape (flat)
Shape cavities	around
Type plate	Microtiter plate
Type frame MTP	without frame

