



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Brand Deep Well Plate 96 Well non-sterile, stackable, 1.2 ml, PP, Low profile

**€179.60**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

Stackable deep well plates are one of the most important components in cell culture with many applications: cultivation of microorganisms, extraction of nucleic acids and proteins, sample storage or screening studies.

- DNA-free, DNase & RNase-free
- temperature range from -80 to 121°C
- resistant to DMSO

## Additional Information

No.	BN-1516
Manufacturer (Brand)	Brand
old neoLab article no.	7-5226
EAN	4058072208738
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	White
Material	Polypropylene (PP)
Volume	1.2 ml
sterile	No
Multi-hard disk features	alphanumeric coded
Number of cavities	96
Fluid behavior	hydrophobic
Shape floor	U-shape (round)
Shape cavities	around
Profile	Low Profile
Type cavity rim	no edge (flat)
Type plate	Deep Well Plates
Type frame MTP	without frame

