



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Brand Deep Well Plate 96 Well non-sterile, stackable, 0.5 ml, PP

€221.70
plus VAT &
Shipping

Product Images



Description

Stackable deep well plates are one of the most important components in cell culture with many applications: cultivation of microorganisms, extraction of nucleic acids and proteins, sample storage or screening studies.

- DNA-free, DNase & RNase-free
- temperature range from -80 to 121°C
- resistant to DMSO

Additional Information

| | |
|--------------------------|--------------------|
| No. | BN-1518 |
| Manufacturer (Brand) | Brand |
| old neolab article no. | 7-5222 |
| EAN | 4058072208752 |
| Transport temperature | Room temperature |
| autoclavable at 121°C | Yes |
| Color | White |
| Material | Polypropylene (PP) |
| Volume | 0.5 ml |
| sterile | No |
| Multi-hard disk features | alphanumeric coded |
| Number of cavities | 96 |
| Fluid behavior | hydrophobic |
| Shape floor | U-shape (round) |
| Shape cavities | around |
| Profile | Standard |
| Type cavity rim | no edge (flat) |
| Type plate | Deep Well Plates |
| Type frame MTP | without frame |

