



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Eppendorf Pipette Tips, Deepwell plate 96/2000µl, Normal Pack, Standard, white, 20 plates

Special Price

€150.00 was

€200.00

**plus VAT &
Shipping**

Product Images



Description

Deepwell Plate 96/2000 µL, wells clear, 2.000 µL, PCR clean, white, 20 plates (5 bags × 4 plates) - OptiTrack® matrix: up to 30% faster sample retrieval and fewer pipetting errors due to high-contrast alphanumeric labeling - RecoverMax® well design: Optimized well geometry for maximum sample recovery and excellent mixing properties - Minimal residual volume and high homogeneity from well to well - Raised well edges and smooth surface for reliable sealing, even during heat sealing - Easy and reliable stacking - g-Safe®: exceptional centrifugation resistance for up to 6.000 × g for faster protocols and improved sample quality - High-quality polypropylene for high resistance to chemicals, mechanical stress and extreme temperatures - Manufactured without demolding aids, plasticizers and biocides - Substances that can negatively affect bioassays - Autoclavable (121 °C, 20 min) - Available with barcode. For further information see: (URL: <https://www.eppendorf.com/barcode>)[www.eppendorf.com/barcode]

Additional Information

No.	VB-0493
Manufacturer (Brand)	Eppendorf
EAN	4058072187156
Transport temperature	Room temperature
Color	White
Volume	2000 µl
Number of cavities	96
Fluid behavior	hydrophobic

