



neoLab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**neoCulture Multitest Plates 24 well, flat bottom,  
PS, sterile, individually sterile packed, box of 100  
pcs**

Special Price

**€107.00** was

**€139.00**

**plus VAT &  
Shipping**

## Product Images



## Description

---

Cell and tissue culture plates with treated surface are perfect for adherent cell growth in multiple, comparative and other analysis. Uniform well volume ensures equal growth surface. The well surface is smooth and free of streaking to maximize usable growth area. To reduce evaporation, raised edges are provided on each well with matching uniform rings on the inside of the lid. The single position (cut corner) lid reduces the risk of cross-contamination and handling hazards. Wells are marked with alphanumeric code for easy identification. The plates are suitable for all common instruments and automats.

Sterilized by gamma radiation

DNase / Rnase free

Non-pyrogenic

## Additional Information

---

No.	C-8208
Manufacturer (Brand)	neoCulture
EAN	4058072154523
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (crystal clear)
Material	Polystyrene (PS)
Volume	3.5 ml
sterile	Yes
Multi-hard disk features	surface optimized individually packed alphanumerically coded
Number of cavities	24
Fluid behavior	hydrophobic
Shape floor	F shape (flat)
Shape cavities	around
Type plate	Cell culture plates
Type frame MTP	full frame
suitable for cooling in	Refrigerator freezer/ freezer

