



neoLab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neoCulture Multitest Plates 12-well flat, PS, sterile, individually packed, box of 100 pieces


Special Price

€107.00 was

€139.00

plus VAT &
Shipping

Product Images



Description

Cell and tissue culture plates with treated surface are perfect for adherent cell growth in multiple, comparative and other analysis. Uniform well volume ensures equal growth surface. The well surface is smooth and free of streaking to maximize usable growth area. To reduce evaporation, raised edges are provided on each well with matching uniform rings on the inside of the lid. The single position (cut corner) lid reduces the risk of cross-contamination and handling hazards. Wells are marked with alphanumeric code for easy identification. The plates are suitable for all common instruments and automats.

Sterilized by gamma radiation

DNase / Rnase free

Non-pyrogenic

Additional Information

No.	C-8209
Manufacturer (Brand)	neoCulture
EAN	4058072154707
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (crystal clear)
Material	Polystyrene (PS)
Volume	6.8 ml
sterile	Yes
Multi-hard disk features	surface optimized individually packed alphanumerically coded
Number of cavities	12
Fluid behavior	hydrophobic
Shape floor	F shape (flat)
Shape cavities	around
Type plate	Cell culture plates
Type frame MTP	full frame
suitable for cooling in	Refrigerator freezer/ freezer

