



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz ASTM centrifuge tube, cylindrical, graduation up to 100 ml

€39.46
plus VAT &
Shipping



Description

ASTM centrifuge glass, cylindrical

With conical lower part, according to ASTM D 91

Capacity: 100 ml

loadable up to RZB 700 (relative centrifuge loading)

Length: max. 203 mm

Additional Information

No.	LZ-1581
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3584
EAN	4058072313197
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	37 mm
Color	transparent (clear)
Material	Borosilicate glass DURAN®
Volume	100 ml
sterile	No
Height	203 mm
Reaction vessels and test tubes Properties	graduated
Floor outside	conical. pointed
Type edge	Flared rim
Type Tube/ Micro Vessel	Centrifuge tubes
Wall thickness	normal

