



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz ASTM centrifuge tube, cylindrical,
graduation up to 100 ml, marks at 50 ml and
100 ml**

**€40.91
plus VAT &
Shipping**

Product Images



Description

ASTM centrifuge tube, cylindrical

With conical lower part, according to ASTM D 96 (old)

Capacity: 100 ml

Loadable up to RZB 700 (relative centrifuge loading)

Length: max. 167 mm

Additional Information

No.	LZ-1582
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3585
EAN	4058072313203
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	44 mm
Color	transparent (clear)
Material	Borosilicate glass DURAN®
Volume	100 ml
sterile	No
Height	167 mm
Reaction vessels and test tubes Properties	graduated
Floor outside	conical. pointed
Type edge	Flared rim
Type Tube/ Micro Vessel	Centrifuge tubes
Wall thickness	normal

