



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e


Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® Autoanalyzer Cup 3.5 ml, f. Olympus, 1000 pcs./pack

**€27.90**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

made of polystyrene, transparent, for the most common autosamplers.  
1000 pieces per pack

## Additional Information

---

No.	E-1272
Manufacturer (Brand)	neoLab
EAN	4058072063764
Transport temperature	Room temperature
autoclavable at 121°C	No
diameter	15 mm
Color	transparent (crystal clear)
Material	Polystyrene (PS)
Volume	3.5 ml
Height	49 mm
TBST MAX	70 °C
TBST MIN	-20 °C
Reaction vessels and test tubes Properties	for 1xuse
Floor outside	conical. with stand ring
Type edge	Flared rim
Type Tube/ Micro Vessel	Test tubes/tubes
Wall thickness	strong-walled

