




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Autoanalyzer Cup 0.25 ml, ES300/60, 1000 pcs./pack

**Price on
request
plus VAT &
Shipping**



Description

Autoanalyzer cups made of polystyrene, transparent, for the most common autosamplers.

Volume: 0.25 ml

1.000 pieces per pack.

Additional Information

No.	E-1273
Manufacturer (Brand)	neoLab
EAN	4058072076696
Transport temperature	Room temperature
autoclavable at 121°C	No
diameter	17 mm
Color	transparent (crystal clear)
Material	Polystyrene (PS)
Volume	3 ml
Height	38 mm
TBST MAX	70 °C
TBST MIN	-20 °C
Reaction vessels and test tubes Properties	for 1xuse
Floor outside	conical. with stand ring
Type edge	Flared rim
Type Tube/ Micro Vessel	Test tubes/tubes
Wall thickness	strong-walled

