



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



DURAN® NMR tubes, for four application frequencies with caps, 178 x 4.95 mm +/- 0.05 mm

€193.00 As low
as
€148.50
plus VAT &
Shipping

Product Images



Description

NMR test tubes are characterized by tightest tolerances and highest precision.
(straightness, wall thickness and wall thickness distribution)

Application: for all standard NMR measurements.

Made of borosilicate glass 3.3 with labeling field and lid made of PP.

Economic: 300 MHz

Additional Information

No.	SD-0939
Manufacturer (Brand)	DURAN
old neoLab article no.	E-3990
EAN	4058072171872
Transport temperature	Room temperature
autoclavable at 121°C	No
diameter	4,95 mm
Color	transparent (clear)
sterile	No
Length	178 mm
Floor outside	around
Type Tube/ Micro Vessel	NMR tubes
Wall thickness	thin-walled

