



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Nalgene PPCO Centrifuge Bottles with Sealing Cap Case of 24 13.700 x G

€151.40 As low
as
€121.84
plus VAT &
Shipping

Product Images

Description

Nalgene™ PPCO Centrifuge Bottles with Seal Cap Centrifuge large volumes of biological and chemical samples at low to medium speed with these PPCO centrifuge bottles equipped with a seal cap. Capacity (metric): 450 mL Capacity (anglo): 16.9 oz. Color: Translucent Material: polypropylene copolymer (PPCO) Number Per Package: 4 Type: centrifuge bottle Nominal Capacity: 500 mL Closure Included: Yes Autoclavable: Yes Closure Material: polypropylene (PP) Closure Type: seal closure unit with silicone gasket Graduated: No Height (Anglo-American): 6.27 in. Max. Force: 13,700 x G Length (Anglo): 2.735 in. Shape: Round Dispensing Through Closure: No Width (Anglo-American): 2.735 in. Unit Size: Case of 24

Additional Information

No.	CD-0976
Manufacturer (Brand)	Nalgene
VGKL number	238731412
EAN	4058072169664
Gross weight	2
Transport temperature	Room temperature
Material	Polypropylene (PP)
Volume	500 ml
Country of origin	USA
Packaging volume	0.03400000 ccm
Customs tariff number	39233010

