



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Nalgene PPCO centrifuge bottles with seal closure Case of 36 27,500 x g

€115.04 As low
as
€92.68
plus VAT &
Shipping

Product Images



Description

Nalgene™ PPCO Centrifuge Bottles with Seal Cap Centrifuge large volumes of biological and chemical samples at low to medium speed with these PPCO centrifuge bottles equipped with a seal cap. Capacity (metric): 250 mL Capacity (anglo): 8.45 oz. Color: Translucent Material: polypropylene copolymer (PPCO) Number Per Package: 4 Type: centrifuge bottle Nominal Capacity: 250 mL Closure Included: Yes Autoclavable: Yes Closure Material: Polypropylene (PP) Closure Type: Silicone Gasket Cap Set Graduated: No Height (Anglo-American): 5.23 in. Max. Force: 27,500 x g Length (anglo): 5.23 in. Shape: Round Dispensing Through Closure: No Width (Anglo-American): 2.40 in. Unit Size: Case of 36

Additional Information

No.	CD-0975
Manufacturer (Brand)	Nalgene
VGKL number	238731415
EAN	4058072172824
Gross weight	2.245
Transport temperature	Room temperature
Material	Polypropylene (PP)
Volume	250 ml
Country of origin	USA
Packaging volume	0.02800000 cc
Customs tariff number	39233010

