




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® Reagent tubes 5 ml, round, 75 x 12 mm,
1000 pcs./pack**

**€28.70
plus VAT &
Shipping**

Product Images



Description

Reagent and centrifuge tubes made of polypropylene.

Ideally suited for a wide range of laboratory work, round bottom, high thermal and mechanical resistance, high resistance to organic and inorganic reagents, autoclavable at 121°C/2.1 bar.

Material:

Polypropylene

Additional Information

No.	E-1649
Manufacturer (Brand)	neolab
EAN	4058072983086
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	12 mm
Color	transparent (clear)
Material	Polypropylene (PP)
Volume	5 ml
sterile	No
Height	75 mm
TBST MAX	121 °C
TBST MIN	-20 °C
Reaction vessels and test tubes Properties	with special chemical resistance
Floor outside	around
Type edge	straight edge
Type Tube/ Micro Vessel	Test tubes/tubes
Wall thickness	normal

