



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e


Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**neolab® Reagent tubes 12 ml, round, 100 x 16 mm, 2000 pcs./pack**

**€114.20  
plus VAT &  
Shipping**

## Product Images



## Description

---

Reagent and centrifuge tubes made of polypropylene.

Ideally suited for a wide range of laboratory work, round bottom, high thermal and mechanical resistance, high resistance to organic and inorganic reagents, autoclavable at 121°C/2.1 bar.

**Material:**

Polypropylene

## Additional Information

---

No.	E-1650
Manufacturer (Brand)	neoLab
EAN	4058072066253
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	16 mm
Color	transparent (clear)
Material	Polypropylene (PP)
Volume	12 ml
sterile	No
Height	100 mm
Reaction vessels and test tubes Properties	with special chemical resistance
Floor outside	around
TBST MAX	121 °C
TBST MIN	-20 °C
Type edge	straight edge
Type Tube/ Micro Vessel	Test tubes/tubes
Wall thickness	normal

