



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Hecht test tubes AR glass, 11.6 mm Ø, 75 mm lg., 200 pcs./pack

Special Price


€12.60 was

€14.43 As low as

€11.41

**plus VAT &
Shipping**

Product Images



Description

AR-Glas® (soda glass)

Original Hecht-ASSISTENT

Strong-walled, approx. 0.7 mm

Round base, smooth rim.

Additional Information

No.	WQ-2396
Manufacturer (Brand)	Hecht
old neoLab article no.	E-2406
EAN	4058072443740
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	11,6 mm
Color	transparent (clear)
Material	Soda-lime glass (normal glass) 1
Volume	5 ml
Height	75 mm
TBST MAX	121 °C
TBST MIN	-196 °C
Reaction vessels and test tubes Properties	with special mechanical strength with special chemical resistance
Floor outside	around
Type edge	straight edge
Type Tube/ Micro Vessel	Test tubes/tubes
Wall thickness	strong-walled

