



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Hecht test tubes AR glass, 14.5 mm Ø, 135 mm lg., 100 pcs./pack

Special Price


**€17.30** was

**€20.15** As low as

**€16.53**

**plus VAT &  
Shipping**

### Product Images



## Description

---

AR-Glas® (soda glass)

Original Hecht-ASSISTENT

Strong-walled, approx. 0.7 mm

Round base, smooth rim.

## Additional Information

---

No.	WQ-2402
Manufacturer (Brand)	Hecht
old neoLab article no.	E-2412
EAN	4058072443757
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	14,5 mm
Color	transparent (clear)
Material	Soda-lime glass (normal glass) 1
Volume	14 ml
Height	135 mm
TBST MAX	121 °C
TBST MIN	-196 °C
Reaction vessels and test tubes Properties	with special mechanical strength with special chemical resistance
Floor outside	around
Type edge	straight edge
Type Tube/ Micro Vessel	Test tubes/tubes
Wall thickness	strong-walled

