



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Hecht test tubes AR glass, 15.25 mm Ø, 100 mm lg., 100 pcs./pack

Special Price


€10.70 was

€12.51 As low as

€9.89

**plus VAT &
Shipping**

Product Images



Description

AR-Glas® (soda glass)

Original Hecht-ASSISTENT

Strong-walled, approx. 0.7 mm

Round base, smooth rim.

Additional Information

No.	WQ-2403
Manufacturer (Brand)	Hecht
old neoLab article no.	E-2413
EAN	4058072443764
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	15,25 mm
Color	transparent (clear)
Material	Soda-lime glass (normal glass) 1
Volume	12 ml
Height	100 mm
TBST MAX	121 °C
TBST MIN	-196 °C
Reaction vessels and test tubes Properties	with special mechanical strength with special chemical resistance
Floor outside	around
Type edge	straight edge
Type Tube/ Micro Vessel	Test tubes/tubes
Wall thickness	strong-walled

