



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Hecht test tubes, thick-walled, with smooth rim No. 42775009 Durham test tube -

€9.55 As low as
€7.56
plus VAT &
Shipping

Product Images

Description

Test tubes, thick-walled, with smooth rim mm Ø approx. 40x 8 No. 42775009 Durham test tube - Tube Durham No. 42775024 Kahn glass, Coombs test tube - Kahn tube, Coombs test tube - Tube de Kahn, Coombs tube à essais No. 42775027 für Blutgruppenuntersuchung - for blood test - pour l'examen de grou Old item number: 2775/2

Additional Information

No.	WQ-2392
Manufacturer (Brand)	Hecht
VGKL number	111814003
EAN	4058072753344
Gross weight	0.164
Transport temperature	Room temperature
diameter	8 mm
Material	Soda-lime glass
Country of origin	Germany
Customs tariff number	70179000
Length	40 mm

