



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz round-bottomed glass with standard ground joint, without graduation, with PE stopper, Ø 18 x H. 150 mm

**€70.70
plus VAT &
Shipping**

Product Images



Description

Test tubes, with standard ground joint, without graduation
Round-bottomed glass with NS sleeve, with PE stopper.

Additional Information

No.	LZ-1603
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3606
EAN	4058072313418
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	18 mm
Color	transparent (clear)
Material	Borosilicate glass DURAN®
Volume	30 ml
sterile	No
Height	150 mm
Floor outside	around
Standard ground sleeve	14/23
Type edge	Standard ground
Type Tube/ Micro Vessel	Test tubes/tubes
Wall thickness	normal

