



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz round-bottomed glass with standard ground joint, without graduation, with PE stopper, Ø 28 x H. 200 mm

€16.08 As low as
€12.37
plus VAT & Shipping

Product Images



Description

Test tubes, with standard ground joint, without graduation
Round-bottomed glass with NS sleeve, with PE stopper.

Additional Information

No.	LZ-1607
Manufacturer (Brand)	Lenz
old neolab article no.	B-3610
EAN	4058072313456
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	28 mm
Color	transparent (clear)
Material	Borosilicate glass DURAN®
Volume	100 ml
sterile	No
Height	200 mm
Floor outside	around
Standard ground sleeve	24/29
Type edge	Standard ground
Type Tube/ Micro Vessel	Test tubes/tubes
Wall thickness	normal

