



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Lenz round bottom glasses with NS 14/23, 100 x  
16 mm, 10 pcs./pack**

**€40.30  
plus VAT &  
Shipping**

## Product Images



## Description

---

DURAN<sup>®</sup>, sleeve 14/23, without stopper,  
10 pieces per pack

## Additional Information

---

No.	LZ-1592
Manufacturer (Brand)	Lenz
old neolab article no.	B-3595
EAN	4058072313302
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	16 mm
Color	transparent (clear)
Material	Borosilicate glass
Height	100 mm
TBST MAX	121 °C
TBST MIN	-196 °C
Reaction vessels and test tubes Properties	with special chemical resistance
Floor outside	around
Standard ground sleeve	14/23
Type edge	Standard ground
Type Tube/ Micro Vessel	Test tubes/tubes
Wall thickness	normal

