



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz round bottom glasses with NS 14/23, 150 x
18 mm, 10 pcs./pack**

€8.15 As low as

€6.27

**plus VAT &
Shipping**

Product Images



Description

DURAN®, sleeve 14/23, without stopper,
10 pieces per pack

Additional Information

No.	LZ-1595
Manufacturer (Brand)	Lenz
old neolab article no.	B-3598
EAN	4058072313333
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	18 mm
Color	transparent (clear)
Material	Borosilicate glass
Height	150 mm
TBST MAX	121 °C
TBST MIN	-196 °C
Reaction vessels and test tubes Properties	with special chemical resistance
Floor outside	around
Standard ground sleeve	14/23
Type edge	Standard ground
Type Tube/ Micro Vessel	Test tubes/tubes
Wall thickness	normal

