



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Brand tubes 1.2 ml, palletized, PP, rack with 96
single tubes, non-sterile**

**€115.85
plus VAT &
Shipping**

Product Images



Description

The autoclavable tubes can be used in many ways in the laboratory, whether for PCR, EIA, storage or transport.

The vials are freezable down to -80°C .

The lids made of PE are not suitable for autoclaving.

Additional Information

No.	BN-2281
Manufacturer (Brand)	Brand
old neoLab article no.	7-5238
EAN	4058072214906
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	9 mm
Color	transparent (crystal clear)
Material	Polypropylene (PP)
Volume	1.2 ml
sterile	No
Height	45 mm
TBST MAX	121°C
TBST MIN	-80°C
Reaction vessels and test tubes Properties	can be used as cryotubes
Floor outside	around
Floor inside	around
suitable for cooling in	Ultra Deep Freezer/ Deep Freezer
Type Tube/ Micro Vessel	Reaction vessels

