



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e


Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Brand tubes 1.2 ml, palletized, PP, tubes single, non-sterile

**€71.75**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

The autoclavable tubes can be used in many ways in the laboratory, whether for PCR, EIA, storage or transport.

The vials are freezable down to  $-80^{\circ}\text{C}$ .

The lids made of PE are not suitable for autoclaving.

## Additional Information

---

No.	BN-2283
Manufacturer (Brand)	Brand
old neoLab article no.	7-5240
EAN	4058072214920
Transport temperature	Room temperature
autoclavable at $121^{\circ}\text{C}$	Yes
diameter	9 mm
Color	transparent (crystal clear)
Material	Polypropylene (PP)
Volume	1.2 ml
sterile	No
Height	45 mm
TBST MAX	$121^{\circ}\text{C}$
TBST MIN	$-80^{\circ}\text{C}$
Reaction vessels and test tubes Properties	can be used as cryotubes
Floor outside	around
Floor inside	around
suitable for cooling in	Ultra Deep Freezer/ Deep Freezer
Type Tube/ Micro Vessel	Reaction vessels

