



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Brand tubes 1.2 ml, palletized, PP, rack with 12
strips of 8 single tubes, non-sterile**

**€115.25
plus VAT &
Shipping**

Product Images



Description

The autoclavable tubes can be used in many ways in the laboratory, whether for PCR, EIA, storage or transport.

The vials are freezable down to -80°C .

The lids made of PE are not suitable for autoclaving.

Dimensions of the tubes:

- Outer diameter: 8.8 mm

- Height: 45 mm

Additional Information

No.	BN-2282
Manufacturer (Brand)	Brand
old neoLab article no.	7-5239
EAN	4058072214913
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	9 mm
Color	transparent (crystal clear)
Material	Polypropylene (PP)
Volume	1.2 ml
sterile	No
Height	45 mm
Reaction vessels and test tubes Properties	can be used as cryotubes
Floor outside	around
Floor inside	around
TBST MAX	121 °C
TBST MIN	-80 °C
Type Tube/ Micro Vessel	Reaction vessels
suitable for cooling in	Ultra Deep Freezer/ Deep Freezer

