



neoLab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Brand tubes 1.2 ml, palletized, PP, strips of 8  
tubes, non-sterile, loose**

**€71.75  
plus VAT &  
Shipping**

## Product Images



## Description

---

The autoclavable tubes can be used in many ways in the laboratory, whether for PCR, EIA, storage or transport.

The vials are freezable down to  $-80^{\circ}\text{C}$ .

The lids made of PE are not suitable for autoclaving.

Dimensions of the tubes:

- Outer diameter: 8.8 mm

- Height: 45 mm

## Additional Information

---

No.	BN-2284
Manufacturer (Brand)	Brand
old neoLab article no.	7-5241
EAN	4058072214937
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	9 mm
Color	transparent (crystal clear)
Material	Polypropylene (PP)
Volume	1.2 ml
sterile	No
Height	45 mm
TBST MAX	121 °C
TBST MIN	-80 °C
Reaction vessels and test tubes Properties	can be used as cryotubes
Floor outside	around
Floor inside	around
suitable for cooling in	Ultra Deep Freezer / Deep Freezer
Type Tube / Micro Vessel	Reaction vessels

